VR3 SEATED LEG CURL TOTAL ACCESS

FEATURES

Five position adjustable calf pad accommodates different lower leg lengths.

Large seat base facilitates easy transfer from a wheelchair.

Thigh pad hold-down mechanism allows the thighs to be secured above the knee for reduced shear force on the knee joint while keeping it in a stable position.

Raised numbers and large font on ROM. Start positions are easy to adjust/identify for visually impaired and older users.

Start ROM on cam synchronized with user.

Regardless of the starting position, the Cybex patent pending ROM adjustment ensures the cam moves to maintain the correct strength profile for the user.

User accessible weight stack offers ease of weight selection from seated exercise position.



SPECIFICATIONS

Product Number	14061
Dimensions L × W × H	56" \times 39" \times 64" (142 cm \times 99 cm \times 162 cm)
Standard Weight Stack	165 lbs (75 kg)
Heavy Weight Stack	205 lbs (93 kg)



Thigh pad clamping mechanism



Raised Numbering



Upgrade your club with customized Cybex equipment. Visit mydesign.cybexintl.com





Large seat base

WORLD HEADQUARTERS

10 Trotter Drive \cdot Medway \cdot MA 02053 USA T + 1.508.533.4300 \cdot F + 1.508.533.5500 cybexintl.com

CYBEX INTERNATIONAL UK LTD Premier House · Beveridge Lane Bardon Hill · Coalville · LE67 1TB UK T + 44.845.606.0228 · F + 44.845.606.0227 CYBEX PRODUCTS ARE DESIGNED AND BUILT IN THE USA





 \odot 2015, Cybex International, all rights reserved. Specifications subject to change 04/15. Actual color may differ from samples shown.